


# Description chart of surgical sutures codes-1



## 1. Common



1 THREAD		
RASTORO	PGA (Polyglycolic Acid)	G
RASTORO FAST	PGA Rapid	K
RADIK	PGLA	A
RADIK FAST	PGLA Rapid	J
GILZA	PDO (Polydioxanone)	D
SIMFRA	PGCL	Q
KOLLGUT PLAIN	Plain Gut	L
KOLLGUT CHROM	Chromic Gut	C
SILDA	Silk	S
MOLYLON	Monofilament Nylon	N
MOLYLON	Braided Nylon	B
Supramid	Supramid	V
PIDILEN	Polypropylene	P
RETEL	Polyester	E
KOFLON	PTFE	Z
KOLLNOVA	PVDF	PV
RIKID	Steel	GS
KOLLFORCE	UHDPE	H

2 USP SIZE		
1	1-0	
2	2-0	
3	3-0	
4	4-0	
5	5-0	
6	6-0	
7	7-0	
8	8-0	
9	9-0	
A	10-0	
X	11-0	
Y	12-0	
B	1	
C	2	
D	3	
E	4	
F	5	

3 NEEDLE TIP		
Taper	1	
Blunt Point	2	
Reverse Cutting	3	
Cutting	4	
Diamond	5	
Premium Cutting	6	
Taper-cut	7	
Spatula/Lancet	8	
Square	9	
Wing	A	
Sleigh Taper	B	
CC	C	
Trocar	T	

4 NEEDLE LENGTH	
Needle below 10mm, needle length can be decimal, such as 6.4mm, unit: mm	

5 NEEDLE BODY		
1/2	2	
3/5	3	
1/4	4	
5/8	5	
Compound	6	
Straight	0	
J	J	


6 ADDITIONAL INFO	
Additional item (needle before thread)	

7 THREAD LENGTH	
Thread length (1-3 numbers) unit: cm	


8 THREADS PER ENVELOPE	
Common package with multi-pieces x product number	

Thread	Default color		
PGA	Violet	W	White
PGLA	Violet	G	Green
PDO	Violet	B	Black
RPGLA	Undyed	V	Violet
RPGA	Undyed	L	Blue
PGCL	Undyed	X	Red
E	Green	U	Undyed
N	Blue (Ophthalmic black)	H	Heavy needle
PTFE	White	P	Loop
P	Blue	R	Reel
S	Black	A	Black
PVDF	Blue	D	Double
H	Mixed color	T	Pledget
V	Undyed	K	Clip
		N	Uncoated

## 2. Precut




## 3. Double needle with one thread




Mark: If same needle shape, should follow short to long.

## 4. Double needle with two thread



## 5. Cassette



Mark: If same needle shape, should according to the needle length