

Description chart of surgical sutures codes-2



I. With Needle

1		2		3		4		5		6		7		8			9		10		
Type of Thread		Thread Color		Size of Thread		Needle shape		Wire Size		Length of needle		Needle Curve		Thread length			Needle quantity		Accessories		
PL	POLYGLACTIN 910 (PGLA)	RADIK	W	White	USP	Metric	TP	Taper Point (Taper Body)	H	Heavy wire size	03	3 mm	A	1/4 circle	inch	inch	metric	0	0	L	LOOP
PF	POLYGLACTIN Rapid	RADIK FAST	B	Black	#	All	TC	Taper Cut (Taper body)	N	Normal wire size	04	4 mm	B	1/2 circle	013	5"	13cm	1	1	EL	TAIL LOOP
PY	POLYGLYCOLIC ACID (PGA)	RASTORO	L	Blue	120	12-0	BP	Blunt Point (Taper body)	T	Slim/Thin Wire size	05	5 mm	C	3/8 circle	020	8"	20cm	2	2	P	Pledget
PR	POLYGLYCOLIC ACID Rapid	RASTORO FAST	G	Green	110	11-0	RC	Reverse Cutting (Triangular body)	A	Visi-black	06	6 mm	D	5/8 circle	030	12"	30cm			B	Barbed
PG	POLYGLYCAPRONE 25	SIMFRA	V	Violet	100	10-0	CC	Conventional Cutting (Square body)	F	Ski Taper	6.5	6.5 mm	E	Straight	045	18"	45cm			R	Reel
PD	POLYDIOXANONE	GILZA	U	Undyed	090	9-0	SP	Conventional Spatula (Flat body)	E	Ski cutting	08	8 mm	G	Bicurve	050	20"	50cm			K	Button
PS	POLYESTER	RESTEL	M	Mixed Color	080	8-0	RP	Reverse Cutting - Prime	G	Ski reverse cutting	09	9 mm	J	Hook	060	24"	60cm			Z	tape
PP	POLYPROPYLENE	PIDILEN	X	Pink	070	7-0	CP	Conventional Cutting Prime	R	ROTATING	10	10 mm			070	27"	70cm				
NM	NYLON POLYAMID Monofilament	MOLYLON	Y	yellow	060	6-0	DD	Diamond			11	11 mm			075	29"	75cm				
NB	NYLON POLYAMID Multifilament	MOLYLON	N	UNCOATED	050	5-0	TR	Trocar			12	12 mm			090	36"	90cm				
SB	SILK BRAIDED	SILDA			040	4-0	SQ	Square			13	13 mm			100	39"	100cm				
GC	CHROMIC GUT	KOLLGUT CHROM			030	3-0					14	14 mm			120	47"	120cm				
GP	PLAIN GUT	KOLLGUT PLAIN			020	2-0				 mm			125	49"	125cm				
SS	STAINLESS STEEL	RIKID			000	0					99	99 mm			135	53"	135cm				
PLA	POLYGLACTIN 910 (PGLA) Antibacterial	RADIK PLUS			001	1					00	100 mm			150	59"	150cm				
PVA	POLYGLYCOLIC ACID Antibacterial	RASTORO PLUS			002	2									180	71"	180cm				
PGN	POLYGLYCAPRONE 25 Antibacterial	SIMFRA PLUS			003	3									200	79"	200cm				
PDA	POLYDIOXANONE Antibacterial	GILZA PLUS			004	4															
PT	PTFE	KOFLON			005	5															
EPT	ePTFE	e-KOFLON			006	6															
PV	PVDF	KOLLNOVA																			
FC	FORCE FIBER	KOLLFORCE																			

II Without needle

1		2		3		4			5		
Type of Thread		Thread Color		Size of Thread		Thread length			Thread quantity		
PL	POLYGLACTIN 910 (PGLA)	RADIK	W	White	USP	Metric	inch	inch	metric	Qty.	
PF	POLYGLACTIN Rapid	RADIK FAST	B	Black	#	All	013	5"	13cm	01	1
PY	POLYGLYCOLIC ACID (PGA)	RASTORO	L	Blue	120	12-0	020	8"	20cm	02	2
PR	POLYGLYCOLIC ACID Rapid	RASTORO FAST	G	Green	110	11-0	030	12"	30cm	03	3
PG	POLYGLYCAPRONE 25	SIMFRA	V	Violet	100	10-0	045	18"	45cm	04	4
PD	POLYDIOXANONE	GILZA	U	Undyed	090	9-0	050	20"	50cm	05	5
PS	POLYESTER	RESTEL	M	Mixed Color	080	8-0	060	24"	60cm	06	6
PP	POLYPROPYLENE	PIDILEN			070	7-0	070	27"	70cm	07	7
NM	NYLON POLYAMID Monofilament	MOLYLON			060	6-0	075	29"	75cm	08	8
NB	NYLON POLYAMID Multifilament	MOLYLON			050	5-0	090	36"	90cm	09	9
SB	SILK BRAIDED	SILDA			040	4-0	100	39"	100cm	10	10
GC	CHROMIC GUT	KOLLGUT CHROM			030	3-0	120	47"	120cm	11	11
GP	PLAIN GUT	KOLLGUT PLAIN			020	2-0	125	49"	125cm	12	12
SS	STAINLESS STEEL	RIKID			000	0	135	53"	135cm	13	13
PLA	POLYGLACTIN 910 (PGLA) Antibacterial	RADIK PLUS			001	1	150	59"	150cm	14	14
PVA	POLYGLYCOLIC ACID Antibacterial	RASTORO PLUS			002	2	180	71"	180cm	15	15
PGN	POLYGLYCAPRONE 25 Antibacterial	SIMFRA PLUS			003	3	200	79"	200cm	16	16
PDA	POLYDIOXANONE Antibacterial	GILZA PLUS			004	4					
PT	PTFE	KOFLON			005	5					
EPT	ePTFE	e-KOFLON			006	6					
PV	PVDF	KOLLNOVA									
FC	FORCE FIBER	KOLLFORCE									

Code sample For KOLLSUT needle Suture:

PLV030TPH36B070-2

Description:

1	Type of thread:	Polyglactin 910	PL
2	Color of thread:	Violet	V
3	Size of thread:	UPS 3-0	030
4	Needle shape:	Taper point	TP
5	Wire size:	Heavy	H
6	Length of needle:	36 mm	36
7	Needle type:	1/2 circle	B
8	Length of thread:	70 cm	070
9	Needle quantity:	2 needles (Loop)	2

Code sample for KOLLSUT non needle suture:

PLV03006013

Description:

1	Type of thread:	Polyglactin 910	PL
2	Color of thread:	Violet	V
3	Size of thread:	UPS 3-0	030
4	Length of thread:	60 cm	060
5	Thread quantity:	13 threads	13

Code sample for KOLLSUT non needle suture:

PS020TCN17B075-2P1X10

Description:

1	Type of thread:	Polyester	PS
2	Size of thread:	UPS 2-0	020
3	Needle shape:	Taper cut	TP
4	Wire size:	Normal	N
5	Length of needle:	17 mm	17
6	Needle type:	1/2 circle	B
7	Length of thread:	75 cm	075
8	Needle quantity:	2 needles (Loop)	2
9	Pledget		P
10	1	3x6x1.5	
	2	3x3x1.5	
	11	Thread quantity: 10 threads	